

Allora - Shortcut Key Reference Guide

The following list gives you an overview of the function key combinations:

 **Stop any sound from playing. (This is helpful when if you made make a mistake.)**

 **The 'message' key is used to recall a stored message.**

Press the 'message' key first and then press a figure or a letter. The message that is associated with this key will then be called. This also applies to the combination:

- message key + Fn + figure or character
- message key + Shift + figure or character
- message key + Fn + Shift + figure or character

To save a message under a key, first type the message, then press the 'message' key  twice followed by the figure or the letter under which you want to save the message. Then press 'confirm'  or 'ignore' . You can also use the following combinations as locations to save a message:



- Fn + figure or character
- Shift + figure or character
- Fn + Shift + figure or character

 **The 'sound' key is used to play stored sounds or recordings.**

Press the 'sound' key first, and then press a figure or a letter to the play the associated recording. This also applies to the combination:

- sound key + Fn + figure or character
- sound key + Shift + figure or character
- sound key + Fn + Shift + figure or character

Pressing the sound key twice will make the 'attention' sound.

	Increase the volume.
	Decrease the volume.
	Go to the Allora settings menu.
	Mute / unmute the sound.
	Add the message to the sentence prediction list.
	Pronounce date.
	Pronounce time.
	Increase the character size. Once maximum size is reached, the next step will start from the minimum size and then continues increasing again.
	Decrease the character size. Once minimum size is reached, the next step will start from the maximum size and then continues decreasing again.
	Switch the 2nd display on or off. This setting is stored until it is changed again.
	Switch between sound from the 2nd display or from the Allora (only when 2nd display is disconnected). This setting is stored until it is changed again.
	Switch between word prediction, sentence prediction and no prediction. This setting is stored until it is changed again.
	Show the Allora battery status, the 2nd display status and the version of the firmware.
	Show the telecom provider, the signal strength of the phone operator (texting only) and the free memory for text messages.

